

ARSA

ENTERPRISES

Ophthalmic Instruments

Speculums

21-3162

Pediatric Speculum

- Spring with solid blade
- Overall length 35mm



21-3163



Berraquer Speculum

- Solid blade, 14mm, angled, temporal
- Overall length 40mm

21-3164

Kratz-Berraquer Speculum

- Wire blade, 15mm, angled, temporal
- Overall length 40mm



21-3165

Berraquer Speculum

- Angled, temporal
- Overall length 40mm, width 30mm

21-3166

Leibinger Speculum

- Solid blade, 15mm
- Angled, temporal
- Overall length 65mm



21-3167

Leibinger Speculum

- Wire blade, 14mm, angled, temporal
- Overall length 65mm



21-3168

Neo-Natal Speculum

- Wire blade, 6mm
- Angled, temporal
- Overall length 65mm




Capsulorrhexis forceps

21-3169

Utrata Capsulorrhexis Forceps


- Round handle, direct action
- 45° Straight shaft
- 0.4mm Sharp tip
- Overall length 103mm



21-3170

Utrata Capsulorrhexis Forceps

- Round handle, direct action
- Curved shaft
- 0.4mm Sharp tip
- Overall length 103mm



21-3171

Inamura Capsulorrhexis Forceps

- Round handle, direct action
- Curved shaft
- 0.4mm Sharp tip
- Overall length 105mm



21-3172

Utrata Capsulorrhexis Forceps


- Flat handle, direct action
- Curved shaft
- 0.4mm Sharp tip
- Overall length 105mm



21-3173

Utrata Capsulorrhexis Forceps


- Flat handle, direct action
- 45° Straight shaft
- 0.4mm Sharp tip
- Overall length 105mm



21-3174

Mics Capsulorrhexis Forceps

- Micro-incision, 1.8mm
- Round handle
- Curved shaft, direct action
- 0.4mm Sharp tip
- Overall length 97mm



21-3175

Utrata Capsulorrhexis Forceps

- Flat handle
- Curved shaft, direct action
- 0.4mm Sharp tip
- Overall length 105mm



21-3176

Fulcher Capsulorrhexis Forceps

- Flat handle
- Curved shaft, direct action
- 0.4mm Rounded tip
- Overall length 105mm



21-3177

Utrata Capsulorrhexis Forceps

- Flat handle
- 45° shaft, direct action
- Overall length 105mm



21-3178

Mics Capsulorrhexis Forceps

- Round handle
- 1.8mm incision, direct action
- 45° Shaft
- Overall length 105mm




Tissue forceps

21-3179

Bonn Tissue Forceps

- 1x2 teeth 0.25mm
- 6mm Tying platform
- Overall length 84mm



21-3180

Bonn Micro Tissue Forceps

- 1x2 teeth 0.25mm
- 6mm Tying platform
- Overall length 106mm



21-3181

ST. Martins Tissue Forceps


- 1x2 teeth 0.6mm
- Overall length 90mm



21-3182

ST. Martins Tissue Forceps


- 1x2 teeth 0.3mm
- Overall length 90mm



21-3183

Lister Utility Forceps

- 1x2 teeth
- 6mm Tying platform
- Overall length 114mm



21-3184

ST. Martins Tissue Forceps


- 0.3mm Notch
- 5mm Tying platform
- Overall length 108mm



21-3185

Mcpherson Tying Forceps


- Straight plain tip
- 6mm Tying platform
- Overall length 85mm



21-3186

Mcpherson Tying Forceps

- Straight plain tip
- 6mm Tying platform
- Overall length 110mm



21-3187

Max Fine Tying Forceps

- Curved smooth tip
- No platform
- Overall length 110mm



21-3188

Kelman-Mcpherson Tying Forceps


- 8mm Tying platform, angled
- Overall length 85mm



21-3189

Colibri Toothed Forceps

- 1x2 teeth 0.25mm
- Overall length 65mm



21-3190

Colibri Toothed Forceps

- 1x2 teeth 0.25mm
- Overall length 110mm




Utility forceps

21-3191

Epilation Forceps


- Angled tips
- Overall length 91mm



21-3192

Moorfields Suture & Utility Forceps


- 14mm Serrated tying platform
- Overall length 107mm



21-3193

Mosquito Forceps

- Ring handle
- Straight tip
- Overall length 100mm




Picture is smaller than actual size

21-3194

Mosquito Forceps

- Ring handle
- Curved tip
- Overall length 100mm



Picture is smaller than actual size



21-3195

Jewelers Forceps

- 45° angle, 7mm tip
- Overall length 110mm



21-3196

Adson Toothed Forceps

- Dissecting 0.25mm teeth
- Overall length 120mm



21-3197

Adson Non-Toothed Forceps

- Overall length 120mm
- Serrated

21-3198

Adson-Brown Forceps

- Overall length 120mm



21-3199

Jewelers No.5 Forceps

- Overall length 110mm



21-3200

Watchmaker No.7 Forceps

- 13mm tip
- Overall length 115mm



Manipulators / Choppers

21-3201

Sinsky Lens Hook

- Straight
- Tip 0.5mm
- Overall length 124mm



21-3202

Lester IOL Manipulator

- Tip angled at 45°
- Angled at 12mm from tip
- 1mm Branch
- Overall length 125mm



21-3203

Y-Shaped Nucleus Manipulator


- On vertical axis
- Tip angled at 45°
- Angled at 13mm from tip
- Overall length 124mm



21-3204

Neuhann Spatula

- Vaulted spatula angled 11mm from tip



21-3205

Ambidextrous Chopper

- Tip angled at 45°
- Angled at 11mm from tip
- Overall length 115mm



21-3206

Sadiq Chopper

- Tip angled at 45°
- Angled at 11mm from tip
- Overall length 121mm



21-3207

Claw Chopper

- Tip angled at 45°
- Angled at 11mm from tip
- Distal tip length 1mm
- Overall length 120mm



21-3208

Mushroom-Tipped Manipulator

- Tip angled at 45°
- Angled at 11mm from tip
- Overall length 123mm



21-3209

**Sinskey/Mushroom-Tip
Double Ended**

- Tip angled at 45°
- Sinskey tip length 0.65mm
- Overall length 123mm



21-3210

**Aycliff/Sadiq Double
Ended Chpper**

- Tip angled at 45°
- Aycliff tip: 0.5mm
- Sadiq tip: 1.5mm
- Overall length 120mm



21-3211

**Mushroom/Sadiq Double
Ended Chopper**

- Tip angled at 45°
- Sadiq tip: 1.5mm
- Overall length 143mm



21-3212

Sinskey Lens Hook

- Tip angled at 45°
- Tip length: 0.5mm
- Overall length 120mm



21-3213

Drysdale Paddle

- Tip angled at 45°
- Angled tip 11mm from tip
- Spatula 1mm
- Overall length 120mm




Needle holders

21-3214

Castroviejo Needle Holder

- Straight with lock
- 12mm Tying platform
- Overall length 135mm



21-3215

Castroviejo Needle Holder

- Straight without lock
- 12mm Tying platform
- Overall length 135mm



21-3216

Castroviejo Needle Holder

- Delicate curve without lock
- 12mm Tying platform
- Overall length 135mm



21-3217

Troutman Needle Holder

- Curved without lock
- 6mm Tying platform
- Overall length 105mm



21-3218
Caliper & Rule



21-3219
Double Ended Marker

- 3.5mm & 4.0mm ends



21-3220
Thornton Ring

- Globe fixation ring
- Multiple blunt teeth



21-3221
Degree Gauge

- 11.7mm inside diameter
- Marked in 5° increments



21-3222
Alignment Marker

- 3, 6, 9 O'clock placement



21-3223
Toric IOL Marker

- 11mm inside diameter



Scissors

21-3224

Westcott Scissors

- Straight tip
- Overall length 120mm



21-3225

Vannas Scissors

- Straight
- Overall length 95mm



21-3226

Iris Scissors

- Straight, Rounded tip
- Overall length 100mm



21-3227

Iris Scissors

- Straight, Sharp tip
- Overall length 100mm



21-3228

Westcott Scissors

- Curved tip, Sharp
- Overall length 120mm



21-3229


Westcott Scissors

- Curved tip, Rounded
- Overall length 120mm

21-3230

Westcott Scissors

- Dolphin nose tip
- Overall length 120mm



Scissors

21-3231

Westcott Scissors Small

- Curved tip, Sharp
- Short blades
- Overall length 120mm



21-3232


ONG Scissors

- Angled, sharp tip
- Overall length 92mm

21-3233

Vannas Scissors

- Curved tip, Sharp
- Overall length 95mm



21-3234


Iris Scissors

- Curved, sharp tip
- Overall length 100mm

21-3235

Strabismus Scissors

- Curved, rounded tip
- Overall length 115mm



21-3236

Strabismus Scissors

- Straight, rounded tip
- Overall length 115mm

21-3237

Dressing Stitch Scissors

- Sharp/Blunt tip
- Overall length 115mm



Muscle hooks & Probes

21-3238

Graefe Strabismus Hook

- Medium, tip length 10mm
- Overall length 140mm



21-3239

Alexandra Muscle Hook

- Tip length 11mm
- Overall length 130mm



21-3240

Liebreich Probe

- 0000 & 000
- Overall length 125mm



21-3241

Liebreich Probe

- 00 & 0
- Overall length 125mm



21-3242

Liebreich Probe

- 1 & 2
- Overall length 125mm



21-3243

Liebreich Probe

- 3 & 4
- Overall length 125mm



21-3244

Chalazion Clamp

- Small
- Clamp diameter 11mm
- Overall length 90mm



21-3245

Chalazion Clamp

- Medium
- Clamp diameter 18mm x 11mm
- Overall length 90mm

21-3246

Chalazion Clamp

- Large
- Clamp diameter 23mm x 15mm
- Overall length 90mm



21-3247

Desmarres Clamp

- Entropion left
- Clamp diameter 24mm x 16mm
- Overall length 95mm

21-3248

Desmarres Clamp

- Entropion right
- Clamp diameter 24mm x 16mm
- Overall length 95mm



21-3249

Bulldog Clamp

- Overall length 53mm

Curettes & Dilators

21-3250

Curette

- Double ended
- 2mm & 3mm cup
- Overall length 130mm



21-3251

Curette

- Double ended
- 2mm & 5mm cup
- Overall length 130mm

21-3252

Curette

- Double ended
- 2mm & 4mm cup
- Overall length 130mm



21-3253

Punctal Dilator

- Double ended
- Overall length 138mm



21-3254

Nettleships Dilator

- Single ended
- Overall length 125mm



21-3255

Epithelial Spatula

- Spatula length 7mm
- Overall length 113mm



21-3256

Alcohol Well

- Inside diameter 9mm
- Overall length 114mm



21-3257

Bard-Parker Scalpel Handle #3


- Overall length 135mm



21-3258

Backhaus Towel Clip

- Overall length 95mm



Picture is smaller than actual size

21-3259

Desmarres Retractor

- 14mm Width
- Overall length 130mm

